

Step 1

Download attachment file and keep into /home folder

Step 2

Open terminal (ctrl+alt+T) and type following command

```
sudo apt-get install unrar
```

```
cd /usr/local
```

```
sudo tar --strip-component 1 -xvf /home/$USER/node-v4.2.2-linux-x64.tar.gz
```

Step 3

Goto /home folder and extract WeighScaleinNode.zip

Step 4

In terminal type following command

```
cd ~/WeighScaleinNode/
```

```
npm install serialport
```

```
sudo chmod 0777 /dev/ttyUSB0
```

Step 5

TO MAKE USER OF WM EVERY TIME OPEN TERMINAL RUN BELOW TWO COMMANDS AND USE LINK:

```
cd WeighScaleinNode/
```

```
node server.js
```